Summary of Safety and Clinical Performance (SSCP)

of

Medical device 1xklysma saline

Information intended for users and healthcare professionals: Please refer to chapter I

Information intended for patients or lay persons: Please refer to chapter II

Table of Content

I	Summary of Safety and Clinical Performance for healthcare professionals	4
1	Device identification and general information.....	4
1.1	Device trade names	4
1.2	Manufacturer's name and address.....	4
1.3	Manufacturer's single registration number (SRN).....	4
1.4	Basic UDI-DI.....	4
1.5	Medical device nomenclature description/text.....	4
1.6	Class of device	5
1.7	Year when the first certificate (CE) was issued covering the device.....	5
1.8	Authorised representative	5
1.9	Notified Body's name and single identification number	5
2	Intended use of the device	5
2.1	Intended purpose	5
2.2	Indication and target population	5
2.3	Contraindication and/or limitations	6
3	Device description.....	6
3.1	Description of the device	6
3.2	Previous generation and its differences	7
3.3	Description of any accessories.....	7
3.4	Description of any other devices and products	7
4	Risks and warnings.....	8
4.1	Residual risks and undesirable effects	8
4.2	Warnings and precautions.....	8
4.3	Other relevant aspects of safety	9
5	Summary of clinical evaluation and post-market clinical follow-up (PMCF)	9
5.1	Summary of clinical data related to equivalence device	9
5.2	Summary of clinical data from conducted investigations of the device before the CE-marking	9
5.3	Summary of clinical data from other sources	10
5.4	An overall summary of clinical performance and safety	10
5.5	Ongoing or planned post-market clinical follow-up	10
6	Possible diagnostic or therapeutic alternatives.....	11
7	Suggested profile and training for users.....	11
8	Applied Harmonised Standards and CS	12
II	Summary of safety and clinical performance for lay persons	13
9	Device identification and general information.....	13
9.1	Device trade name.....	13
9.2	Manufacturer's name and address.....	13
9.3	Basic UDI-DI.....	13

9.4	Year when the device was first CE-marked.....	13
10	Intended use of the device	14
10.1	Intended purpose	14
10.2	Indications and intended patient groups.....	14
10.3	Contraindications	14
11	Device description.....	15
11.1	Device description and material/substances in contact with patient tissues	15
11.2	Information about medicinal substances in the device	15
11.3	Description of how the device is achieving its intended mode of action.....	15
11.4	Description of accessories.....	15
12	Risks and warnings.....	15
12.1	How potential risks have been controlled or managed	15
12.2	Remaining risks and undesirable effects.....	16
12.3	Warnings and precautions.....	16
12.4	Summary of any field safety corrective action, (FSCA including FSN)	17
13	Summary of clinical evaluation and post-market clinical follow-up	17
13.1	Clinical background of the device	17
13.2	The clinical evidence for the CE-marking	17
13.3	Safety	18
14	Possible diagnostic or therapeutic alternatives.....	18
14.1	General description of therapeutic alternatives.....	18
15	Suggested training for users	18
	Revision history	19

I Summary of Safety and Clinical Performance for healthcare professionals

This Summary of Safety and Clinical Performance (SSCP) is intended to provide public access to an updated summary of the main aspects of the safety and clinical performance of the device.

The SSCP is not intended to replace the Instructions For Use as the main document to ensure the safe use of the device, nor is it intended to provide diagnostic or therapeutic suggestions to intended users or patients.

The following information is intended for users and healthcare professionals.

Following this information there is a summary intended for patients.

1 Device identification and general information

1.1 Device trade names

The trade names of the medical device Klysma are:

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1.2 Manufacturer's name and address

Manufacturer according to MPDG and MDR 2017/745:

Serumwerk Bernburg AG

Hallesche Landstraße 105 b

D- 06406 Bernburg

Germany

1.3 Manufacturer's single registration number (SRN)

DE-MF-000005169

1.4 Basic UDI-DI

Basic Unique Device Identification (Basic UDI-DI): 426020097KLYSMA8M

1.5 Medical device nomenclature description/text

EMDN code: G99 (GASTROINTESTINAL DEVICES – OTHER)

1.6 Class of device

Classification according Annex VIII of Regulation (EU) 2017/745, rule 21: III

1.7 Year when the first certificate (CE) was issued covering the device

The product was first classified as class I medical device acc. to rule 5 MDD in year 2017. Thus, the conformity assessment procedure has been performed by the manufacturer acc. to Annex VII of MDD. As no assessment of the CE conformity by a Notified Body was required, a CE certificate covering the device was not issued in the past.

1.8 Authorised representative

not applicable

1.9 Notified Body's name and single identification number

mdc medical device certification GmbH

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70191 Stuttgart

Germany

notified body number: 0483

2 Intended use of the device

2.1 Intended purpose

For rapid stool evacuation of the rectum, e.g. before operations, diagnostic procedures such as X-ray examinations, colonoscopies or childbirth.

Please refer to the instruction for use for further information.

2.2 Indication and target population

For fast emptying of stool from the lower part of the intestine, the rectum, e.g. before surgery, diagnosis procedures such as x-ray examinations or colonoscopies and births.

It is used without any other device or medicinal products. These products are used either by health care professionals or by patients.

2.3 Contraindication and/or limitations

Contraindications:

- Patients with pre-existing electrolyte disorders
- Patients with impaired kidney function
- Patients with cardiac insufficiency
- Intestinal movement disorder due to an obstruction or blockage
- Conditions in which the intestines are inflamed or swollen and that can lead to cramps, pain, constipation and diarrhoea, e.g. in Hirschsprung's disease, ulcerative colitis, megacolon
- Blood in the stool, stomach pain, abdominal pain
- Nausea, vomiting
- Injury to the anal region caused, for example, by a perforation of the bowel or a fissure
- Postoperative inability to evacuate stool through the bowel
- Conditions in which it is likely that increased levels of electrolytes (salts in the blood), such as sodium or phosphate, are absorbed from the intestine
- Hyperphosphataemia
- Simultaneous intake of preparations containing sodium phosphate (also orally)
- Allergies to the ingredients
- Infants and children under 6 years of age

3 Device description

3.1 Description of the device

1xklysma saline is a rectal sodium phosphate containing enema, filled in PE container (135 ml) with a special shape, suitable for the intended use (rectal administration).

The solution contained in the device is not sterile and incorporates sodium benzoate as preservative.

The shelf life of 1xklysma saline is 36 months. The product should be stored at normal storage conditions.

One or ten containers are packed in a printed carton with instruction for use insert.

1xklysma saline is intended for single use only.

The composition of the solution is the following (for 100 ml):

Sodium dihydrogen phosphate 2 H ₂ O	16.00 g
Disodium hydrogen phosphate 2 H ₂ O	2.98 g
Sodium benzoate	0.12 g
Purified water	to 100 ml

3.2 Previous generation and its differences

The medical device is unchanged in its form and composition.

3.3 Description of any accessories

There are no applicable accessories for the medical device according to Regulation (EU) 2017/745 Article 2 (2). Please note, that you have to lightly grease the tip of enema (e.g. with Vaseline).

3.4 Description of any other devices and products

There are no other devices and products which are intended to be used in combination with the medical device.

4 Risks and warnings

4.1 Residual risks and undesirable effects

Residual risk “Overdose”

An overdose can cause increased blood phosphate and decreased blood calcium or other serious electrolyte imbalances that can lead to dehydration, cardiovascular problems or kidney failure.

Adverse reactions

1xklysma is well tolerated when used as indicated. However, adverse reactions possibly associated with 1xklysma saline have been infrequently reported. In rare cases, adverse events may occur, especially if the medical device is misused.

- Irritation of the intestinal mucosa can occur.
- In rare cases, loss of soft tissue around the rectum (colonic gangrene) due to severe ulceration or wounds has been associated with the use of phosphate-containing enemas in elderly patients.
- It can cause a high level of phosphate or a low level of calcium or other electrolyte imbalances. Dehydration and muscle cramps can result from the water and electrolyte imbalances.

All serious incidents and adverse reactions related to the product must be reported to the manufacturer and to the competent authority of the member state where the user and/or patient is established.

4.2 Warnings and precautions

Precaution:

Use with caution in:

- Elderly, infirm or frail patients
- Patients suffering from dehydration
- Patients suffering from uncontrolled high blood pressure

Interactions with other agents:

In patients who take diuretics or other medications that affect the electrolyte balance, the likelihood of electrolyte imbalances (high phosphate and sodium levels / low calcium levels) is increased when using phosphate-containing enemas or oral laxatives.

No other sodium phosphate preparations, including oral sodium phosphate solutions or tablets, should be administered at the same time.

Warnings:

- Only intended for single use. The reuse or reprocessing of single-use products can lead to contamination.
- Dispose of plastic containers after use, even if not the whole content has been used. Only administer clear solutions in intact containers.
- Do not use after the expiry date stated on the bottle and folding box.
- Keep out of the reach of children!

4.3 Other relevant aspects of safety

No information concerning serious incidents was received.

No field safety corrective actions were taken.

5 Summary of clinical evaluation and post-market clinical follow-up (PMCF)

5.1 Summary of clinical data related to equivalence device

The equivalence of the comparative products* was demonstrated regarding clinical, biological and technical characteristics. The clinical evaluation based on literature of sodium phosphate containing enemas. The equivalent device's SSCP is not available in Eudamed. The phosphate-containing solutions have a similar composition and are intended for colorectal cleansing. All materials are well-known and are widely used in the medical sector. The comparative devices are used under similar conditions of use, e.g. population and clinical condition, and have similar relevant performance according to the expected clinical effect for the specific intended use. The principle of action of phosphate containing enemas based on osmotic forces. The differences, presented in the Clinical Evaluation Report, are minor and not expected to affect the clinical performance and safety. Clinical data on sodium phosphate enemas as well as decades of experience with the clinical use of 1xklyisma saline allow to state that there are no doubts about the efficacious and safe application of such enemas. Based on the review of scientific literature of sodium phosphate containing enemas it can be concluded that the benefits of 1xklyisma saline outweigh their potential risks.

**Klistier- Fresenius, Fleet Enema- C.B. Fleet Company, Clyssie- B.Braun, Cleen Ready to use- Casen Recortadi*

5.2 Summary of clinical data from conducted investigations of the device before the CE-marking

not applicable

5.3 Summary of clinical data from other sources

Sodium phosphate enemas are products widely used in both in-patient and out-patient settings. The most common indication is for symptomatic treatment of constipation, and to a lesser extent in preparation for colonoscopy or surgery. These products have been widely used for many years, and have been associated with minimal adverse effects in the general population. Water and electrolyte disturbances are the reason for the main side effects occurring in patients administered sodium phosphate enemas. However, It should not be forgotten that sodium phosphate enemas are widely used and there are a minimum number of patients with side effects.

5.4 An overall summary of clinical performance and safety

Summarizing all facts discussed in the clinical evaluation report (CER) 1xklyisma saline have been effectively applied as rectal cleansing for decades.

Based on all critical considerations concerning 1xklyisma saline and comparable devices they are without doubt effective in the claimed indications of application. An assessment of the substances contained in 1xklyisma saline did not reveal any hazardous effects when administered according to the precautions specified in the product information. No adverse events are to be expected using the 1xklyisma saline following the recommendations specified in the instruction for use. Data on similar medical devices as well as decades of experience with the clinical use of 1xklyisma saline allow to state that there are no doubts about the efficacious and safe application of such enemas.

An overall review of the findings allows the conclusion that the instruction for use of 1xklyisma saline adequately reflects the prevailing scientific knowledge on its intended purposes. Under the conditions specified in the instruction for use 1xklyisma saline can be recommended as enema for all applications listed in the product information.

It can be concluded that the benefits of 1xklyisma saline outweigh their potential risks. Potential risks related to the administration of 1xklyisma saline are well known and can generally be avoided by following the manufacturer's instruction for use. The overall residual risks and the benefit/risk ratio of the device are fully acceptable. Therefore, the requirements of the Medical Devices Regulation 2017/745 regarding safety and clinical performance of 1xklyisma saline are fulfilled.

5.5 Ongoing or planned post-market clinical follow-up

1xklyisma saline is a well analysed and clinically proven medical device. The safety and performance of phosphate-containing enemas is demonstrated for decades. All clinically relevant information can be supported by clinical data in the Clinical Evaluation Report. The therapeutic benefit of 1xklyisma saline for

rectal cleansing outweighs the remaining residuals risks. The overall risk with the use of 1xklyisma saline in the intended indication is low.

The benefit/risk profile is positive and therefore the requirement on acceptable benefit/risk profile is fulfilled. Therefore, no Post-Market-Clinical Follow-up (PMCF) is necessary.

6 Possible diagnostic or therapeutic alternatives

Depending on the patient's individual preferences and medical history, various therapeutic alternatives are available to treat constipation. The following alternatives are primarily oral medications.

- Osmotic cathartic
- Nonabsorbed osmotic agent
- Stimulant laxative (contact irritant)
- Combination osmotic/laxative
- Combination nonabsorbed osmotic/laxative

7 Suggested profile and training for users

The medical device can be used by patients as well as health care professionals. No special training is required for the use of the medical device.

8 Applied Harmonised Standards and CS

Applied harmonised standards as well as other specific standards related to 1xklyisma are listed below. Common specifications are currently not available.

Standard/Guideline	Title /Description	Applied fully/partly
General		
DIN EN ISO 9001:2015-11	Quality management systems - Requirements	fully
DIN EN ISO 13485: 2021-12	Medical devices - Quality management systems - Requirements for regulatory purposes	fully
DIN EN ISO 14971:2022-04	Medical devices - Application of risk management to medical devices	fully
ISO/TR 24971:2020-06	Medical devices - Guidance on the application of ISO 14971	fully
ISO/TR 20416:2020-07	Medical devices — Post-Market surveillance for manufacturers	fully
MPDG	Medical Device Law Implementation Act	fully
MDR	Medical Device Regulation (EU) 2017/745	fully
Ph.Eur.	European Pharmacopoeia current version	fully
Biological and clinical evaluation		
DIN EN ISO 10993-1:2021-05	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management system	fully
MEDDEV 2.7/1 rev. 4	Clinical Evaluation: A Guide for Manufacturers and Notified Bodies under Directives 93/42/EEC and 90/385/EEC	fully
MDCG 2019-9 Rev.1 March 2022	Summary of safety and clinical performance A guide for manufacturers and notified bodies	fully
MDCG 2020-5 April 2020	Clinical Evaluation - Equivalence A guide for manufacturers and notified bodies	fully
Labelling		
DIN EN ISO 15223-1:2022-02	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements	fully
DIN EN ISO 20417:2022-03	Information to be supplied by the manufacturer of medical devices	fully
MDCG 2018-1 April 2021	Guidance on BASIC UDI-DI and changes to UDI-DI April 2021	fully
Usability		
DIN EN 62366-1:2021-08 Appendix C	Application of usability engineering to medical devices	fully

II Summary of safety and clinical performance for lay persons

This Summary of Safety and Clinical Performance (SSCP) is intended to provide public access to an updated summary of the main aspects of the safety and clinical performance of the device. The information presented below is intended for patients or lay persons. A more extensive summary of its safety and clinical performance prepared for healthcare professionals is found in the first part of this document.

The SSCP is not intended to give general advice on the treatment of a medical condition. Please contact your healthcare professional in case you have questions about your medical condition or about the use of the device in your situation. This SSCP is not intended to replace an Implant card or the Instructions For Use to provide information on the safe use of the device.

9 Device identification and general information

9.1 Device trade name

The trade names of the medical device Klysma are:

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9.3 Basic UDI-DI

426020097KLYSMA8M

9.4 Year when the device was first CE-marked

The product was first classified as class I medical device acc. to rule 5 of Medical Device Directive (MDD) in year 2017.

10 Intended use of the device

10.1 Intended purpose

For rapid stool evacuation of the rectum, e.g. before operations, diagnostic procedures such as X-ray examinations, colonoscopies or childbirth.

Please refer to the instruction for use for further information.

10.2 Indications and intended patient groups

For fast emptying of stool from the lower part of the intestine, the rectum, e.g. before surgery, diagnosis procedures such as x-ray examinations or colonoscopies and births. It is used without any other device or medicinal products. These products are used either by health care professionals or by patients.

10.3 Contraindications

Always use this medical device exactly as described in the IFU or as your doctor or pharmacist has told you. Check with your doctor or pharmacist if you are not sure in any cases [e.g. therapeutic indications, contraindications, precaution].

Do not use this product in:

- Patients with pre-existing electrolyte disorders
- Patients with impaired kidney function
- Patients with cardiac insufficiency
- Intestinal movement disorder due to an obstruction or blockage
- Conditions in which the intestines are inflamed or swollen and that can lead to cramps, pain, constipation and diarrhoea, e.g. in Hirschsprung's disease, ulcerative colitis, megacolon
- Blood in the stool, stomach pain, abdominal pain
- Nausea, vomiting
- Injury to the anal region caused, for example, by a perforation of the bowel or a fissure
- Postoperative inability to evacuate stool through the bowel
- Conditions in which it is likely that increased levels of electrolytes (salts in the blood), such as sodium or phosphate, are absorbed from the intestine
- Hyperphosphataemia
- Simultaneous intake of preparations containing sodium phosphate (also orally)
- Allergies to the ingredients

- Infants and children under 6 years of age

11 Device description

11.1 Device description and material/substances in contact with patient tissues

1xklysma saline is a rectal phosphate solution for bowel evacuation. The product is applied rectally by a health care professional or by the patient. The solution contains sodium phosphate, which has an osmotic activity to rapidly softens and dilutes the stool.

11.2 Information about medicinal substances in the device

Not applicable

11.3 Description of how the device is achieving its intended mode of action

The influx of the solution increases the water content and volume of the stools, making them softer and easier to pass, and increases the pressure which leads to a spontaneous and fast bowel evacuation.

11.4 Description of accessories

There are no applicable accessories for the medical device. Please note, that you have to lightly grease the tip of enema (e.g. with Vaseline).

12 Risks and warnings

Contact your healthcare professional if you believe that you are experiencing side-effects related to the device or its use or if you are concerned about risks. This document is not intended to replace a consultation with your healthcare professional if needed.

12.1 How potential risks have been controlled or managed

The 1xklysma has been designed in accordance with documented processes to ensure that it is designed, manufactured, packaged, and labelled in accordance with the current state of the art and meets all requirements of the appropriate regulations. All risks identified during these activities were mitigated as far as possible and are considered acceptable in regards to the clinical benefit of the device. Continued review of all Post Market Surveillance is performed to identify any additional risks that may be identified after the device was placed on the market.

12.2 Remaining risks and undesirable effects

1xklyisma is well tolerated when used as indicated. However, adverse reactions possibly associated with 1xklyisma saline have been infrequently reported. In rare cases, adverse events may occur, especially if the medical device is misused.

- Irritation of the intestinal mucosa may occur.
- In rare cases, loss of soft tissue around the rectum (colonic gangrene) due to severe ulceration or wounds has been associated with the use of phosphate-containing enemas in elderly patients.
- It can cause a high level of phosphate or a low level of calcium or other electrolyte imbalances. Dehydration and muscle cramps can result from the water and electrolyte imbalances.

All serious incidents and side effects related to the product must be reported a doctor or pharmacist, the manufacturer and to the competent authority of the member state where the user and/or patient is established.

12.3 Warnings and precautions

Always use this medical device exactly as described in the IFU or as your doctor or pharmacist has told you. Check with your doctor or pharmacist if you are not sure in any cases [e.g. therapeutic indications, contraindications, precaution].

Use with caution in:

- Elderly, infirm or frail patients
- Patients suffering from dehydration
- Patients suffering from uncontrolled high blood pressure

Warnings:

- Only intended for single use. The reuse or reprocessing of single-use products can lead to contamination.
- Dispose of plastic containers after use, even if not the whole content has been used. Only administer clear solutions in intact containers.
- Do not use after the expiry date stated on the bottle and folding box.
- Keep out of the reach of children!

12.4 Summary of any field safety corrective action, (FSCA including FSN)

Not applicable

13 Summary of clinical evaluation and post-market clinical follow-up

13.1 Clinical background of the device

Phosphate-containing enemas have been successfully used in different medical fields for decades for several indications. Sufficient clinical data as well as decades of experience with the clinical use of 1xklyisma saline allow to state that there are no doubts about the efficacious and safe application of such enemas. 1xklyisma saline is considered as a safe medical device with a positive benefit-risk ratio. Due to the long term use, safety and the low risk to patient further clinical investigations are not necessary.

13.2 The clinical evidence for the CE-marking

Phosphate-containing enemas have been successfully used in different medical fields for decades for several indications. The clinical evaluation based on literature of sodium phosphate containing enemas. The equivalence of the comparative products* was demonstrated regarding clinical, biological and technical characteristics. The equivalent device's SSCP is not available in Eudamed. As part of the clinical evaluation, the equivalence of the same / similar products was demonstrated. All products are enemas that are composed of sodium and phosphate salts. The phosphate-containing solutions have a similar composition and are intended for colorectal cleansing. All materials are well-known and are widely used in the medical sector. The comparative devices are used under similar conditions of use, e.g. population and clinical condition, and have similar relevant performance according to the expected clinical effect for the specific intended use. The principle of action of phosphate containing enemas is based on osmotic forces. The differences, showing in the Clinical Evaluation Report, are minor and not expected to affect the clinical performance and safety. Clinical data on sodium phosphate enemas as well as decades of experience with the clinical use of 1xklyisma saline allow to state that there are no doubts about the efficacious and safe application of such enemas. Based on the review of scientific literature of sodium phosphate containing enemas it can be concluded that the benefits of 1xklyisma saline outweigh their potential risks. The general safety and performance requirements of the Medical Device Regulation (EU) 2017/745 are fulfilled.

*Klistier- Fresenius, Fleet Enema- C.B. Fleet Company, Clyssie- B.Braun, Cleen Ready to use- Casen Recortadi

13.3 Safety

The therapeutic benefit of 1xklyisma saline for rectal cleansing outweighs the remaining residuals risks. The overall risk with the use of 1xklyisma saline in the intended indication is low.

Mainly, side-effects described in the literature are due to overdose or application in a contraindicated patient population (children below 6) or patients at risk (elderly), and thus, can be avoided by following the recommendation given in the instruction for use. No adverse events are to be expected with proper administration following the recommendations given by the manufacturer in the instructions for use.

The benefit/risk profile is positive and therefore the requirement on acceptable benefit/risk profile is fulfilled.

14 Possible diagnostic or therapeutic alternatives

When considering alternative treatments, it is recommended to contact your healthcare professional who can take into account your individual situation.

14.1 General description of therapeutic alternatives

When considering alternative treatments, it is recommended that the patient contact the responsible healthcare professional who can take into account the individual situation.

There are diverse alternative treatment methods:

- Osmotic cathartic
- Nonabsorbed osmotic agent
- Stimulant laxative (contact irritant)
- Combination osmotic/laxative
- Combination nonabsorbed osmotic/laxative

15 Suggested training for users

The medical device can be used by patients as well as health care professionals.

No special training is required for the use of the medical device.

Revision history

SSCP revision number	Date issued	Change description	Revision validated by Notified Body
00	2021-08-25	First edition	<input type="checkbox"/> Yes, Validation language: English <input checked="" type="checkbox"/> No, not yet validated by the NB
01	2022-05-12	Update: Addition of information for lay person	<input checked="" type="checkbox"/> Yes, Validation language: English <input type="checkbox"/> No
02	2024-05-17	<p>The SSCP was primarily changed editorially and formally. The changes to the SSCP do not result from a change to the CER, PSUR; PMCF report, RM report or others.</p> <ul style="list-style-type: none"> - Overall: Adaption acc. MDCG 2019-9 Rev1 - Ch. 1: Addition of corresponding chapters for healthcare prof and patients/lay persons - Ch. 1.1 and 9.1: Addition of the Slovenian and Russian name - Ch. 1.7 and 9.4: Correction of the year - Ch. 2.2 and 10.2: Adaption acc IFU Rev 03 - Ch. 2.3 and 10.3: Adaption acc IFU Rev 03 - Ch. 3.1: Shortening the chapter to avoid redundant information - Ch. 3.3 and 11.4: Addition to grease the tip of enema - Ch. 4 and 12: Correction acc IFU Rev 03 - Ch. 6: Removal of guidelines - Ch. 7 and 15: Addition, that no special training is necessary - Ch. 8: Adaption acc. 2.8.A1-Klyisma-Standards rev05 - Ch. 9.2: Addition of address 	<input type="checkbox"/> Yes, Validation language: English <input checked="" type="checkbox"/> No, not required, as only editorial changes

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Norma/směrnice	Název/popis	Uplatnění zcela/částečně
Obecné		
DIN EN ISO 9001:2015-11	Quality management systems - Requirements	zcela
DIN EN ISO 13485: 2021-12	Medical devices - Quality management systems - Requirements for regulatory purposes	zcela
DIN EN ISO 14971:2022-04	Medical devices - Application of risk management to medical devices	zcela
ISO/TR 24971:2020-06	Medical devices - Guidance on the application of ISO 14971	zcela
ISO/TR 20416:2020-07	Medical devices — Post-Market surveillance for manufacturers	zcela
MPDG	Medical Device Law Implementation Act	zcela
MDR	Medical Device Regulation (EU) 2017/745	zcela
Ph.Eur.	European Pharmacopoeia current version	zcela
Biologické a klinické hodnocení		
DIN EN ISO 10993-1:2021-05	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management system	zcela
MEDDEV 2.7/1 rev. 4	Clinical Evaluation: A Guide for Manufacturers and Notified Bodies under Directives 93/42/EEC and 90/385/EEC	zcela
MDCG 2019-9 Rev.1 Březen 2022	Summary of safety and clinical performance A guide for manufacturers and notified bodies	zcela
MDCG 2020-5 Duben 2020	Clinical Evaluation - Equivalence A guide for manufacturers and notified bodies	zcela
Označení na obalu		
DIN EN ISO 15223-1:2022-02	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements	zcela
DIN EN ISO 20417:2022-03	Information to be supplied by the manufacturer of medical devices	zcela
MDCG 2018-1 Duben 2021	Guidance on BASIC UDI-DI and changes to UDI-DI April 2021	zcela
Použitelnost		
DIN EN 62366-1:2021-08 Příloha C	Application of usability engineering to medical devices	zcela

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Standard/smernica	Naziv/opis	Uporaba v celoti/delno
Splošno		
DIN EN ISO 9001:2015-11	Quality management systems - Requirements	v celoti
DIN EN ISO 13485: 2021-12	Medical devices - Quality management systems - Requirements for regulatory purposes	v celoti
DIN EN ISO 14971:2022-04	Medical devices - Application of risk management to medical devices	v celoti
ISO/TR 24971:2020-06	Medical devices - Guidance on the application of ISO 14971	v celoti
ISO/TR 20416:2020-07	Medical devices — Post-Market surveillance for manufacturers	v celoti
MPDG	Medical Device Law Implementation Act	v celoti
MDR	Medical Device Regulation (EU) 2017/745	v celoti
Ph.Eur.	European Pharmacopoeia current version	v celoti
Biološka in klinična ovrednotenja		
DIN EN ISO 10993-1:2021-05	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management system	v celoti
MEDDEV 2.7/1 rev. 4	Clinical Evaluation: A Guide for Manufacturers and Notified Bodies under Directives 93/42/EEC and 90/385/EEC	v celoti
MDCG 2019-9 Rev.1 Marec 2022	Summary of safety and clinical performance A guide for manufacturers and notified bodies	v celoti
MDCG 2020-5 April 2020	Clinical Evaluation - Equivalence A guide for manufacturers and notified bodies	v celoti
Označevanje		
DIN EN ISO 15223-1:2022-02	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements	v celoti
DIN EN ISO 20417:2022-03	Information to be supplied by the manufacturer of medical devices	v celoti
MDCG 2018-1 April 2021	Guidance on BASIC UDI-DI and changes to UDI-DI April 2021	v celoti
Uporabnost		
DIN EN 62366-1:2021-08 Priloga C	Application of usability engineering to medical devices	v celoti

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