Hohlfaserhämodialysator (VitaPES® P HF Serie)

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Hollow Fiber Dialyzer (VitaPES® P HF Serie)

Indications: VitaPES Dialyzer can be used for the hemodialysis treatment of acute and chronic renal failure. When selecting this product for a treatment, the total extracorporeal blood volume (in combination with the dialyzer, the bloodline set and any other accessories) shall not exceed 10% of the patient's blood volume. Therefore the device is not intended for paediatric use.

Contraindications: No absolute contraindication for hemodialysis treatment. Strict monitoring shall be achieved for the patients who have tendencies of blood bleeding or clotting during the treatment. In case of any complication that would affect the stable condition, the treatment shall be terminated. Patients who experience hypersensitivity reactions against one of the materials in use should not be treated with this product.

Side effects: Adverse hypersensitivity reactions or thrombocytopenia can occur in rare cases. If severe complications occur, the hemodialysis treatment must be discontinued and appropriate action taken. In this case, the blood from the extracorporeal system should not be returned to the patient.

Warning:

- The dialysate and the blood should flow in counter-current.
- The maximum flow rate of dialysate is 800 mL/min.
- During the treatment, the transmembrane pressure shall not exceed 500 mmHg.
- The blood flow rate shall not be less than 150 mL/min, but not more than 500 mL/min.
- The preparation for hemodialysis and treatment should be carried out under aseptic operation to avoid infection.
- The Dialyzer shall be used under certain medical supervision. The Dialyzer shall be used only in combination with a hemodialysis machine, extracorporeal circuit (blood lines) and fluids for hemodialysis that all meet the respective international and local standards.
- The blood pathway is sterile and non-pyrogenic. The device is sterilized using irradiation. DO NOT use the expired product. DO NOT use the device if package is damaged or if protective end caps are not in place.
- This dialyzer is for single use and reuse is strictly prohibited. Reprocessing of this product may lead to adverse patient reactions and/or device failure. It should be discarded according to laws and regulations relevant to disposal of infectious medical waste so as to prevent infection.
- If abnormal conditions arise during the dialysis, such as bubbles, foreign body, blood leak, clotting, or hemolysis, proper measures shall be taken according to doctor's advice.
- Color identification:
Arterial: side with red label marking,
Venous: side with blue label marking

Instructions for Installation:

- Take out the dialyzer from the pouch, and check whether the dialyzer and its components are in good condition.
- Insert the dialyzer vertically into the holder, the arterial end downwards for filling / pre-flushing from the arterial side.
- Make sure the dialyzer stays firmly in the holder.

Instructions for Use:

1. Priming

- To ensure safe and effective treatment, the following instructions for preparing the dialyzer should also be followed when using dialysis machines that allow automated filling and rinsing.
- Prepare no less than 500 ml sterile, pyrogen-free physiological rinsing solution (e.g. NaCl 0.9% or online rinsing solution) and add appropriate amount of heparin if needed under the doctor's advice.
- Place the arterial line and the venous line onto the dialysis machine according to the instructions for use of the blood line tubing set.

MD Medical Device

- Connect the arterial line, the venous line and the dialyzer.
- Set the flow rate of the blood pump to 80 to 100 mL/min, rinse the blood lines and the dialyzer thoroughly with the rinsing solution; the solution should flow from bottom to top to minimise the formation of air bubbles during filling. Filling in the opposite direction should be avoided.
- Turn the flow rate up to 200–300 mL/min, connect the dialysate and the dialyzer. Turn the dialyzer so that the venous end is downward, and remove air from the dialysate fluid compartment thoroughly.

Note: Continue pre-rinsing until the blood pathway is deaerated.

2. Anticoagulation

- Carry out the anticoagulation process in accordance with the doctor's instructions.

3. Termination of Treatment

- For reinfusion and completion of the treatment, follow the instruction on the dialysis machine in use.
- Prepare 500 mL rinsing solution and turn off the blood pump.
- Withdraw the arterial line from the patient and connect it to the rinsing solution then set the blood flow rate at 100 mL/min, meanwhile prevent air from the blood line. Intermittently clamp the venous line to return the blood to the patient as much as possible.
- Turn off the blood pump and disconnect the patient's venous line.

Note: Do not turn off the air monitor system before blood returns completely to prevent the air flowing into the patient from the blood line.

Product Performance: This dialyzer has reliable performance, which can be used for hemodialysis. The basic parameters of product performance and the laboratory data of this series will be provided as follows for reference.

Note:

1. The laboratory data of this dialyzer was measured according to the standard ISO 8637-1.
2. These data represent typical in vitro performance. Actual in vivo performance will differ due to the patient's blood composition and clinical settings.

Special Storage Conditions and Methods: Please avoid crash and exposure to rain, snow, and direct sunlight during transportation. Please store it in a well-ventilated indoor place with storage temperature of 0°C–40°C, with relative humidity no more than 80% and without corrosive gas. DO NOT store it in a warehouse together with chemicals and humid articles.

Shelf life: 3 years after the date of manufacture

After-sale service: if any product inquiry is needed, due to manufacturing quality problem, please keep the primary package of the product for product tracing.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the national competent authority.

Main Structure:

Component	Housing	Dialysis membrane	Cover	Sealant	Protective cap
Material	Polypropylene	VitaSulfone® PES membrane	Polypropylene	Polyurethane	Polypropylene

Model	VitaPES			P 18 HF			P 19 HF			P 20 HF			P 22 HF		
Product REF No. / Order No.	7738			7739			7740			7742					
Test Conditions: Q _D = 500 mL/min, temperature: 37°C ± 1°C, Q _F = 10 mL/min															
Clearance / Q _B (mL/min)	200	300	400	200	300	400	200	300	400	200	300	400			
Urea (mL/min)	195	270	316	197	278	325	197	282	335	198	289	345			
Creatinine (mL/min)	191	250	290	196	260	300	197	270	311	198	279	322			
Phosphate (mL/min)	188	255	315	192	269	334	195	281	340	196	288	352			
Vitamin B ₁₂ (mL/min)	146	165	190	156	185	200	167	200	212	172	211	222			
Inulin (mL/min)	119	131	132	130	140	144	142	150	154	154	159	165			
Pressure drop of blood compartment (mmHg)	<60	<95	<100	<50	<85	<90	<45	<75	<80	<40	<70	<75			
Ultrafiltration Coefficient (mL/hr/mmHg) Q _B = 300 mL/min, TMP = 100 mmHg	55			60			66			72					
Priming volume (mL)	100			109			120			138					
Effective membrane area (m²)	1.6			1.8			2.0			2.2					
Sieving coefficients (S.C.) Q _B = 200 mL/min Q _F = 30 mL/min	Inulin: 1.0 ± 0.1 β2-microglobulin: ≥ 0.7 Myoglobin: 0.55 ± 0.1 Albumin: ≤ 0.01														

	Do not re-use		sterile blood path, sterilized using irradiation		Blood path is non-pyrogenic
	Date of manufacture		Manufacturer		Temperature limit
	Batch code		Use-by date		This end up
	Catalogue number		Keep dry		Fragile, handle with care
	Do not use if package is damaged		Keep away from sunlight		Consult instructions for use
	Caution		Humidity limitation		Unique Device Identification

Keep this Instruction for use after all of the products in this carton are used up.

Manufacturer
MTP Medical Technologies GmbH
Kunstseidenstraße 4, 01796 Pirna, Germany
info@mtp-pirna.de

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