VitaPES® HF High Flux Dialysatoren

ALLGEMEINE INFORMATION

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


## VitaPES® HF High Flux Dialyzers

### GENERAL INFORMATION

**Indications:** The VitaPES® HF dialyzer is designed for single use in chronic hemodialysis (HD) and hemodiafiltration (HDF) treatment for adult patients.

**Contraindications:** There are no general contraindications for the use of this dialyzer in hemodialysis or hemodiafiltration. Patients who experience hypersensitivity reactions against one of the materials in use should not be treated with this product (please refer to the table Technical Data for material specifications).

**Side-effects:** Adverse hypersensitivity reactions can occur in rare cases. If severe complications occur, the hemodialysis treatment must be discontinued and appropriate action taken. In this case, the blood from the extracorporeal system should not be returned to the patient.

**Sterility:** The VitaPES® dialyzer is a product with sterile and pyrogen-free blood path. The product is sterilized by radiation (e-beam).

**Storage:** Store product up to and not longer than the expiry date displayed on the product. Store between 0°C and +30°C.

**Anticoagulation:** The administration and dosage of anticoagulants is the responsibility of the physician. Insufficient anticoagulation may cause clotting or obstruction in the capillaries of the dialyzer.

### WARNING

Do not use the product after the expiry date. The expiry date applies only to a correctly stored product in undamaged packaging.

Sterility of the product is guaranteed only when the packaging is undamaged.

The blood inflow and outflow connectors must remain closed until prior to connecting the dialyzer to the blood line connectors.

Preparation of the extracorporeal system and the dialyzer must take place under aseptic conditions. To ensure safe and effective treatment, the following instructions for preparing the dialyzer should also be followed when using dialysis machines that allow automated filling and rinsing.

This dialyzer membrane is highly permeable to water and so it must be used exclusively with volume-controlled devices with regulated equalization. It is essential to comply with the device manufacturers safety instructions when using highly permeable dialyzer membranes. Because of the technical specification of high flux membranes, it is possible for dialysate to enter the blood circulation at low ultrafiltration rates when highly permeable membranes are used. For this reason, the dialysate must always be ultra-pure.

Reprocessing of the dialyzer for multiple use can result in modifications of product components with a possible negative effect on performance and safety of the dialyzer.

### PRIMING/RINSING OF THE BLOOD COMPARTMENT

Priming the blood compartment of the filter must be carried out before the dialysate fluid compartment is filled. Priming direction should be bottom up. The dialyzer is positioned vertically in the holder. Remove the protective caps from the connectors.

Connect the arterial and venous blood line to the filter using aseptic technique. The arterial connection must be placed downwards and the venous connection upwards. Connect the arterial connection to at least 500 ml of sterile, pyrogen-free, physiological solution (e.g. NaCl 0,9 % or online rinsing solution). Connect the venous blood line to a waste container and prime the entire system.

Recommended blood flow rate is approximately 100 ml/min. It is recommended to discard the first 100 ml of the priming solution. Prime with priming solution, ensuring that all the air is removed (to remove air, gently tap venous end of the dialyzer). Ensure levels set correctly in the arterial and venous chambers (recommendation: fill the arterial chamber up to 50 % and the venous chamber up to 80 %). It is advisable to fill the venous chamber shortly towards the end of priming.

Stop the blood pump and now also connect the venous end of the blood line to the rinsing solution and start the blood pump for the circulation. Set pump speed at least 100 ml/min and recirculate for several minutes.

### FILLING THE DIALYSATE FLUID COMPARTMENT

Turn the dialyzer so that the arterial blood inflow is now directed upward and venous side downward. The filter should be maintained in this position during the entire treatment. Connect the dialysate fluid couplings so that the inflow of dialysate fluid is at the bottom and the outflow for dialysate fluid is at the top. Blood and dialysate fluid must flow in a counter-current direction. Now fill the dialysate fluid compartment so that it is completely free of air (dialysate flow rate should be 500 ml/min). Flush the dialysate fluid compartment for several minutes simultaneously with recirculation of the blood side with approximately 100 ml/min. An ultrafiltration rate of 250 ml/h should thereby be set (note: flushing of the dialysate compartment as well as the ultrafiltration can be initiated during the priming of the blood side.)

### HEMODIAFILTRATION

Ensure that substitution lines are primed and connected only with sterile, pyrogen-free and physiological rinsing solution.

### CONNECTION TO THE PATIENT

Connect the arterial bloodline free of air to the patient's arterial access. The venous bloodline should be connected to the waste bag. Now fill the extracorporeal system (blood flow rate approx. 100 ml/min). Discard as much of the priming fluid as possible and stop the blood pump. Connect the venous bloodline to the patient's venous access. In patients with a tendency to hypotension, the reservoir of the extracorporeal system can be replaced with fresh, sterile infusion solution. Treatment as prescribed by the physician can now be started.

### IN THE EVENT OF DIALYZER BLOOD LEAK

Activate bypass and reduce the blood flow rate. Test the dialysate fluid for the presence of blood. If there is evidence of blood, exchange the dialyzer and the extracorporeal system. Do not return the blood. If necessary, administer replacement fluid to the patient to compensate for the blood loss.

### RECOMMENDATION FOR POSTDILUTION HDF

When increasing volume of water is removed, the fluid proportion of the blood is also drastically reduced, e.g. with excessive substitution or ultrafiltration flow. As a result there is an increased risk of thrombus formation, resulting in obstruction of the dialyzer capillaries. In addition, the diffusion characteristics of the dialyzer membrane would also be altered. In general, the maximum substitution flow should not exceed 25 % of the blood flow.

### DISCONTINUING TREATMENT – RETURNING THE PATIENT'S BLOOD

Stop the blood and ultrafiltration pump. Clamp off the arterial bloodline and the arterial access. Disconnect the arterial bloodline and connect it to a sterile, pyrogen-free and physiological solution (e.g. NaCl 0,9 %). Unclamp the arterial bloodline and rinse-back the patient's blood at the desired flow rate. When the required amount of fluid has been returned (approx. 250 ml) stop the pump. Clamp the venous bloodline and clamp the patient's venous access. Discard the dialyzer, bloodlines and fluid lines into an appropriate container immediately after use. Disposal must comply with the applicable national laws.

### CAUTION

In clinical use, varying results may be obtained resulting from different methods of measurement as well as possible variations in the batches of product.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the national competent authority.

**Presentation:** Bag containing one dialyzer.

| TECHNICAL DATA<br>VitaPES® High Flux Dialyzer                                                                                           |                                  | 150 HF             |     | 170 HF          |     | 190 HF          |     | 210 HF          |     |
|-----------------------------------------------------------------------------------------------------------------------------------------|----------------------------------|--------------------|-----|-----------------|-----|-----------------|-----|-----------------|-----|
| REF-No. / Order-No.                                                                                                                     |                                  | 70104215 / 7215    |     | 70108217 / 7217 |     | 70105219 / 7219 |     | 70106221 / 7221 |     |
| In vitro performance und physical data according to ISO 8637-1 (KUF: human blood)                                                       |                                  |                    |     |                 |     |                 |     |                 |     |
| Clearances (ml/min)*                                                                                                                    | Q <sub>8</sub>                   | 200                | 300 | 200             | 300 | 200             | 300 | 200             | 300 |
|                                                                                                                                         | Urea (ml/min)                    | 196                | 271 | 197             | 276 | 197             | 279 | 198             | 282 |
|                                                                                                                                         | Creatinine (ml/min)              | 188                | 246 | 191             | 254 | 193             | 261 | 195             | 268 |
|                                                                                                                                         | Phosphate (ml/min)               | 182                | 232 | 184             | 238 | 186             | 242 | 189             | 249 |
|                                                                                                                                         | Vitamin B <sub>12</sub> (ml/min) | 147                | 173 | 151             | 185 | 154             | 193 | 159             | 198 |
|                                                                                                                                         | Inulin (ml/min)                  | 93                 | 105 | 103             | 115 | 110             | 121 | 113             | 123 |
| Sieving coefficient                                                                                                                     | Inulin                           | 1<br>0.8<br>< 0.01 |     |                 |     |                 |     |                 |     |
| Q <sub>8</sub> = 300 ml/min                                                                                                             | β2-Microglobulin                 |                    |     |                 |     |                 |     |                 |     |
| Q <sub>f</sub> = 60 ml/min                                                                                                              | Albumin                          |                    |     |                 |     |                 |     |                 |     |
| Membrane surface area                                                                                                                   | (m²)                             | 1.5                |     | 1.7             |     | 1.9             |     | 2.1             |     |
| UF-Coefficient                                                                                                                          | (ml/h/mmHg)                      | 60                 |     | 67              |     | 73              |     | 80              |     |
| Priming volume                                                                                                                          | (ml; blood side)                 | 89                 |     | 99              |     | 112             |     | 123             |     |
| Flow and pressure limits                                                                                                                |                                  |                    |     |                 |     |                 |     |                 |     |
| Max. transmembrane pressure                                                                                                             | (mmHg)                           | 500                |     | 500             |     | 500             |     | 500             |     |
| Max. blood flow                                                                                                                         | (ml/min)                         | 500                |     | 500             |     | 500             |     | 500             |     |
| Max. blood side pressure                                                                                                                | (mmHg)                           | 500                |     | 500             |     | 500             |     | 500             |     |
| Membrane material: Purema® H, internal diameter 200 µm, wall thickness 30 µm                                                            |                                  |                    |     |                 |     |                 |     |                 |     |
| Housing: polycarbonate, Potting compound: polyurethane, O-Ring: silicone, blood and dialysate fluid connections according to ISO 8637-1 |                                  |                    |     |                 |     |                 |     |                 |     |

\*Q<sub>b</sub> = blood flow

Q<sub>o</sub> = dialysate flow = 500 ml/min

Q<sub>f</sub> = ultrafiltration flow = 0 ml/min

Purema® is a registered trademark of 3M.

do not use if package is damaged

Medical Device

Unique Device Identifier

0483

Manufacturer  
MTP Medical Technologies GmbH  
Kunstseidenstraße 4  
01796 Pirna  
Germany  
info@mtp-pirna.de

Storage temperature between 0°C and 30°C

Sterile blood path, Sterilisation by radiation

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